

AssayMax™ Canine Haptoglobin ELISA Kit

Assaypro LLC 3400 Harry S Truman Blvd St. Charles, MO 63301 T (636) 447-9175 F (636) 395-7419 www.assaypro.com

For any questions regarding troubleshooting or performing the assay, please contact our support team at support@assaypro.com.

Thank you for choosing Assaypro.

Assay Summary

Step 1. Add 50 μ l of Standard or Sample per well. Incubate 2 hours.

Step 2. Wash, then add 50 μl of Biotinylated Antibody per well. Incubate 1 hour.

Step 3. Wash, then add 50 μl of SP Conjugate per well. Incubate 30 minutes.

Step 4. Wash, then add 50 μ l of Chromogen Substrate per well. Incubate 20 minutes.

Step 5. Add 50 μ l of Stop Solution per well. Read at 450 nm immediately.

Symbol Key

Consult instructions for use.

Assay Template

12								
11								
10								
6								
∞								
7								
9								
ß								
4								
æ								
2								
1								
	A	B	C	D	Е	Ŀ	Ð	т

AssayMax[™] Canine Haptoglobin ELISA Kit

Catalog No. ECH2003-1 Sample insert for reference use only

Introduction

Haptoglobin (HP, zonulin) is a plasma protein with hemoglobin-binding capacity and a plasma glycoprotein that forms a stable complex with hemoglobin to aid the recycling of heme iron. It is a well-known marker of hemolysis (1). High haptoglobin level in plasma was associated with an increased cardiovascular risk in obese men (2), inflammation (3), atherosclerosis (4), and systemic sclerosis (5).

Principle of the Assay

The AssayMax[™] Canine Haptoglobin ELISA (Enzyme-Linked Immunosorbent Assay) Kit is designed for detection of haptoglobin in canine **plasma and serum samples**. This assay employs a quantitative **sandwich enzyme immunoassay** technique that measures canine haptoglobin in approximately 4 hours. A polyclonal antibody specific for canine haptoglobin has been precoated onto a 96-well microplate with removable strips. Haptoglobin in standards and samples is sandwiched by the immobilized antibody and a biotinylated polyclonal antibody specific for canine haptoglobin, which is recognized by a streptavidin-peroxidase (SP) conjugate. All unbound material is washed away and a peroxidase enzyme substrate is added. The color development is stopped and the intensity of the color is measured.

Caution and Warning

- This product is for **Research Use Only** and is not intended for use in diagnostic procedures.
- Prepare all reagents (diluent buffer, wash buffer, standard, biotinylated antibody, and SP conjugate), as instructed, prior to running the assay.
- Prepare all samples prior to running the assay. The dilution factors for the samples are suggested in this insert. However, the user should determine the optimal dilution factor.
- Spin down the SP conjugate vial and the biotinylated antibody vial before opening and using contents.
- The Stop Solution is an acidic solution.
- The kit should not be used beyond the expiration date.

Reagents

- **Canine Haptoglobin Microplate:** A 96-well polystyrene microplate (12 strips of 8 wells) coated with a polyclonal antibody against canine haptoglobin.
- **Sealing Tapes:** Each kit contains 3 precut, pressure sensitive sealing tapes that can be cut to fit the format of the individual assay.
- **Canine Haptoglobin Standard:** Canine haptoglobin in a buffered protein base (80 ng, lyophilized).
- **Biotinylated Canine Haptoglobin Antibody (50x):** A 50-fold concentrated biotinylated polyclonal antibody against canine haptoglobin (120 µl).
- MIX Diluent Concentrate (10x): A 10-fold concentrated buffered protein base (30 ml).
- Wash Buffer Concentrate (20x): A 20-fold concentrated buffered surfactant (30 ml, 2 bottles).
- SP Conjugate (100x): A 100-fold concentrate (80 µl).
- **Chromogen Substrate (1x):** A stabilized peroxidase chromogen substrate tetramethylbenzidine (7 ml).
- Stop Solution (1x): A 0.5 N hydrochloric acid solution to stop the chromogen substrate reaction (11 ml).

Storage Condition

- Upon arrival, immediately store components of the kit at recommended temperatures up to the expiration date.
- Store SP Conjugate and Biotinylated Antibody at -20°C.
- Store Microplate, Diluent Concentrate (10x), Wash Buffer, Stop Solution, and Chromogen Substrate at 2-8°C.
- Unused microplate wells may be returned to the foil pouch with the desiccant packs and resealed. May be stored for up to 30 days in a vacuum desiccator.
- Store Standard at 2-8°C before reconstituting with Diluent and at -20°C after reconstituting with Diluent.

Other Supplies Required

- Microplate reader capable of measuring absorbance at 450 nm
- Pipettes (1-20 μl, 20-200 μl, 200-1000 μl, and multiple channel)
- Deionized or distilled reagent grade water

Sample Collection, Preparation, and Storage

- **Plasma:** Collect plasma using one-tenth volume of 0.1 M sodium citrate as an anticoagulant. Centrifuge samples at 3000 x g for 10 minutes and collect plasma. An 800000-fold sample dilution is suggested into MIX Diluent; however, user should determine optimal dilution factor depending on application needs. The undiluted samples can be stored at -20°C or below for up to 3 months. Avoid repeated freeze-thaw cycles.
- Serum: Samples should be collected into a serum separator tube. After clot formation, centrifuge samples at 3000 x g for 10 minutes and remove serum. An 800000-fold sample dilution is suggested into MIX Diluent; however, user should determine optimal dilution factor depending on application needs. The undiluted samples can be stored at -20°C or below for up to 3 months. Avoid repeated freeze-thaw cycles.

Applicable samples may also include biofluids, cell culture, and tissue homogenates. If necessary, user should determine optimal dilution factor depending on application needs.

	Guidelines for Dilutions of 100-fold or Greater (for reference only; please follow the insert for specific dilution suggested)				
100x			10000x		
A)	4 μl sample : 396 μl buffer (100x) = 100-fold dilution Assuming the needed volume is less than or equal to 400 μl.	A) B)	4 μl sample : 396 μl buffer (100x) 4 μl of A : 396 μl buffer (100x) = 10000-fold dilution Assuming the needed volume is less than or equal to 400 μl.		
1000x			100000x		
A) B)	4 μl sample : 396 μl buffer (100x) 24 μl of A : 216 μl buffer (10x) = 1000-fold dilution	A) B) C)	4 μl sample : 396 μl buffer (100x) 4 μl of A : 396 μl buffer (100x) 24 μl of B : 216 μl buffer (10x) = 100000-fold dilution		
	Assuming the needed volume is less than or equal to 240 μl.		Assuming the needed volume is less than or equal to 240 $\mu l.$		

Refer to Dilution Guidelines for further instruction.

Reagent Preparation

- Freshly dilute all reagents and bring all reagents to room temperature before use.
- MIX Diluent Concentrate (10x): Dilute the MIX Diluent Concentrate 10fold with reagent grade water to produce a 1x solution. When diluting the concentrate, make sure to rinse the bottle thoroughly to extract any

precipitates left in the bottle. Mix the 1x solution gently until the crystals have completely dissolved. Store for up to 30 days at 2-8°C.

 Canine Haptoglobin Standard: Reconstitute the Canine Haptoglobin Standard (80 ng) with 4 ml of MIX Diluent to generate a 20 ng/ml standard stock solution. Allow the vial to sit for 10 minutes with gentle agitation prior to making dilutions. Prepare duplicate or triplicate standard points by serially diluting from the standard stock solution (20 ng/ml) 2-fold with equal volume of MIX Diluent to produce 10, 5, 2.5, 1.25, 0.625, and 0.313 ng/ml solutions. MIX Diluent serves as the zero standard (0 ng/ml). Any remaining stock solution should be stored at -20°C and used within 30 days. Avoid repeated freeze-thaw cycles.

Standard Point	Dilution	[Haptoglobin] (ng/ml)
P1	1 part Standard (20 ng/ml)	20
P2	1 part P1 + 1 part MIX Diluent	10
P3	1 part P2 + 1 part MIX Diluent	5.0
P4	1 part P3 + 1 part MIX Diluent	2.5
P5	1 part P4 + 1 part MIX Diluent	1.25
P6	1 part P5 + 1 part MIX Diluent	0.625
P7	1 part P6 + 1 part MIX Diluent	0.313
P8	MIX Diluent	0.0

- **Biotinylated Canine Haptoglobin Antibody (50x):** Spin down the antibody briefly and dilute the desired amount of the antibody 50-fold with MIX Diluent to produce a 1x solution. The undiluted antibody should be stored at -20°C.
- Wash Buffer Concentrate (20x): Dilute the Wash Buffer Concentrate 20fold with reagent grade water to produce a 1x solution. When diluting the concentrate, make sure to rinse the bottle thoroughly to extract any precipitates left in the bottle. Mix the 1x solution gently until the crystals have completely dissolved.
- SP Conjugate (100x): Spin down the SP Conjugate briefly and dilute the desired amount of the conjugate 100-fold with MIX Diluent to produce a 1x solution. The undiluted conjugate should be stored at -20°C.

Assay Procedure

- Prepare all reagents, standard solutions, and samples as instructed. Bring all reagents to room temperature before use. The assay is performed at room temperature (20-25°C).
- Remove excess microplate strips from the plate frame and return them immediately to the foil pouch with desiccants inside. Reseal the pouch

securely to minimize exposure to water vapor and store in a vacuum desiccator.

- Add 50 μl of Canine Haptoglobin Standard or sample to each well. Gently tap plate to thoroughly coat the wells. Break any bubbles that may have formed. Cover wells with a sealing tape and incubate for 2 hours. Start the timer after the last addition.
- Wash the microplate manually or automatically using a microplate washer. Invert the plate and decant the contents; hit 4-5 times on absorbent material to completely remove the liquid. If washing manually, wash five times with 200 µl of Wash Buffer per well. Invert the plate each time and decant the contents; hit 4-5 times on absorbent material to completely remove the liquid. If using a microplate washer, wash six times with 300 µl of Wash Buffer per well; invert the plate and hit 4-5 times on absorbent material to completely remove the liquid.
- Add 50 µl of Biotinylated Canine Haptoglobin Antibody to each well. Gently tap plate to thoroughly coat the wells. Break any bubbles that may have formed. Cover wells with a sealing tape and incubate for 1 hour.
- Wash the microplate as described above.
- Add 50 µl of SP Conjugate to each well. Gently tap plate to thoroughly coat the wells. Break any bubbles that may have formed. Cover wells with a sealing tape and incubate for 30 minutes. Turn on the microplate reader and set up the program in advance.
- Wash the microplate as described above.
- Add 50 μl of Chromogen Substrate to each well. Gently tap plate to thoroughly coat the wells. Break any bubbles that may have formed. Incubate in ambient light for 20 minutes or until the optimal blue color density develops.
- Add 50 μl of Stop Solution to each well. The color will change from blue to yellow. Gently tap plate to ensure thorough mixing. Break any bubbles that may have formed.
- Read the absorbance on a microplate reader at a wavelength of 450 nm immediately. If wavelength correction is available, subtract readings at 570 nm from those at 450 nm to correct optical imperfections. Otherwise, read the plate at 450 nm only. Please note that some unstable black particles may be generated at high concentration points after stopping the reaction for about 10 minutes, which will reduce the readings.

Data Analysis

- Calculate the mean value of the duplicate or triplicate readings for each standard and sample.
- To generate a standard curve, plot the graph using the standard concentrations on the x-axis and the corresponding mean 450 nm

absorbance (OD) on the y-axis. The best fit line can be determined by regression analysis using log-log or four-parameter logistic curve fit.

• Determine the unknown sample concentration from the Standard Curve and multiply the value by the dilution factor.

Typical Data

• The typical data is provided for reference only. Individual laboratory means may vary from the values listed. Variations between laboratories may be caused by technique differences.

Standard Point	ng/ml	OD	Average OD
P1	20	2.594	2 570
FI	20	2.546	2.370
P2	10	1.806	2.570 1.795 1.045 0.577 0.302 0.170 0.100
FΖ	10	1.784	
Р3	5.0	1.025	1.0/15
FJ	5.0	1.065	1.045
P4	2.5	0.563	0 577
F4	2.5	0.591	0.377
Р5	1.25 0.294 0.310	0.294	0 302
FJ		0.302	
P6	0.625	0.169	0 170
FO		0.171	0.170
Р7	0.313	0.092	0 100
F 7	F7 0.313		0.100
P8	0.0	0.022	0.025
10 0.0		0.028	0.025
Sample: Poo	oled Normal	0.365	0 277
Sodium Citrate P	lasma (800000x)	0.389	0.377

Standard Curve

• The curve is provided for illustration only. A standard curve should be generated each time the assay is performed.

Reference Value

 Normal canine haptoglobin plasma and serum levels range from 0 – 3 mg/ml.

Performance Characteristics

- The minimum detectable dose of canine haptoglobin as calculated by 2SD from the mean of a zero standard was established to be 79 pg/ml.
- Intra-assay precision was determined by testing three plasma samples twenty times in one assay.
- Inter-assay precision was determined by testing three plasma samples in twenty assays.

	Intra-Assay Precision			Inter-Assay Precision		
Sample	1	2	3	1	2	3
n	20	20	20	20	20	20
CV (%)	3.5%	4.5%	3.1%	6.7%	8.4%	7.8%
Average CV (%)	3.7%				7.6%	

Recovery

Standard Added Value	1 – 10 ng/ml	
Recovery %	87 - 114%	
Average Recovery %	97%	

Linearity

• Plasma samples were serially diluted to test for linearity.

Average Percentage of Expected Value (%)				
Sample Dilution	Plasma			
200000x	106%			
800000x	97%			
320000x	97%			

Cross-Reactivity

Species	Cross-Reactivity (%)
Human	None
Bovine	None
Monkey	None
Mouse	None
Rat	None
Swine	None
Rabbit	None

Troubleshooting

Issue	Causes	Course of Action		
	Use of improper components	 Check the expiration date listed before use. Do not interchange components from different lots. 		
_	Improper wash step	 Check that the correct wash buffer is being used. Check that all wells are empty after aspiration. Check that the microplate washer is dispensing properly. If washing by pipette, check for proper pipetting technique. 		
cisio	Splashing of reagents while loading wells	• Pipette properly in a controlled and careful manner.		
Low Precision	Inconsistent volumes loaded into wells	 Pipette properly in a controlled and careful manner. Check pipette calibration. Check pipette for proper performance. 		
Г	Insufficient mixing of reagent dilutions	 Thoroughly agitate the lyophilized components after reconstitution. Thoroughly mix dilutions. 		
	Improperly sealed microplate	 Check the microplate pouch for proper sealing. Check that the microplate pouch has no punctures. Check that three desiccants are inside the microplate pouch prior to sealing. 		

_	Microplate was left unattended between	 Each step of the procedure should be performed uninterrupted. 			
gna	steps	uninterrupted.			
Unexpectedly Low or High Signal Intensity	Omission of step	 Consult the provided procedure for complete list of steps. 			
or Hi	Steps performed in incorrect order	Consult the provided procedure for the correct order.			
ly Low ol Intensity	Insufficient amount of	 Check pipette calibration. 			
ter Lo	reagents added to wells	 Check pipette for proper performance. 			
⊒ ⊆	Wash step was skipped	 Consult the provided procedure for all wash steps. 			
tec	Improper wash buffer	 Check that the correct wash buffer is being used. 			
xpect	Improper reagent preparation	 Consult reagent preparation section for the correct dilutions of all reagents. 			
Une	Insufficient or prolonged incubation periods	 Consult the provided procedure for correct incubation time. 			
Deficient Standard Curve Fit	Non-optimal sample dilution	 Sandwich ELISA: If samples generate OD values higher than the highest standard point (P1), dilute samples further and repeat the assay. Competitive ELISA: If samples generate OD values lower than the highest standard point (P1), dilute samples further and repeat the assay. User should determine the optimal dilution factor for samples. 			
ndard	Contamination of reagents	 A new tip must be used for each addition of different samples or reagents during the assay procedure. 			
ent Stai	Contents of wells evaporate	 Verify that the sealing film is firmly in place before placing the assay in the incubator or at room temperature. 			
Defici	Improper pipetting	 Pipette properly in a controlled and careful manner. Check pipette calibration. Check pipette for proper performance. 			
	Insufficient mixing of reagent dilutions	 Thoroughly agitate the lyophilized components after reconstitution. Thoroughly mix dilutions. 			

References

- (1) Van Vlierberghe H et al. (2004) Clin Chim Acta. 345(1-2):35-42.
- (2) Engstrom G et al. (2004) Arterioscler Thromb Vasc Biol. 24(8):1498-502.
- (3) Rocha-Pereira P et al. (2004) Br J Dermatol. 150(5):917-28.
- (4) Matuszek MA et al. (2003) Atherosclerosis. 168(2):389-96.
- (5) Kucharz EJ et al. (2000) Clin Rheumatol. 19(2):165-6.

Version 1.9R1